



BEESWIFT
FOCUSED ON SAFETY

Product Code:

RA94640


Selling UOM:

Each

Pack Size:

1

 Direct To Skin

 Fold-to-Pop Mechanism (FTP)

 Instant Activation

 One Time Use

 Warm

SPECIFICATION SHEET

INSTANT WARM PERINEAL COMPRESS 5"X 16"



The Rapid Relief Instant Warm Perineal Compress (Direct-to-skin) quickly delivers soothing warmth, helping to encourage relaxation of tissues and increased blood circulation in the perineal area prior to childbirth, providing maximum patient comfort with minimal nursing intervention.

Features

- Temperature: Delivers approximately 45.6°C/114°F of warm therapy when activated at room temperature
- Extra absorbent compress comes with an adhesive strip to keep it firmly in place
- Longer front-to-back area for better coverage and protection from leakage
- Much more convenient than a warm sitz bath
- FTP (fold-to-pop) technology, allowing activation of the warm compress simply by folding it in half
- Direct-to-skin application

Application

- INSTANT WARM THERAPY RELIEF:
- Encourage relaxation of tissues prior to childbirth
- Start the healing process after childbirth and vasectomy procedures (usually 72 hours after the use of cold therapy)
- Muscle tension and menstrual cramps



Packaging Waste	Type	Units	Weight(g)	Length(mm)	Width(mm)	Height(mm)
Cardboard	Outer	1	260.00	338.00	208.00	312.00
Plastic	Inner	24	4.39			
Carton Weight:	7.20kg					